

# *Unlock* a Path Forward for Your Patients with Urge Urinary Incontinence

# *Revi*<sup>™</sup>

*could be the key  
to relief for patients  
looking for a new  
solution.*



# Demonstrated Efficacy & Safety

## The OASIS Pivotal Trial

A prospective, single-arm, open-label study with 151 patients at 23 medical centers globally was conducted to evaluate the safety and efficacy of the Revi™ System for treating patients with Urge Urinary Incontinence (UUI).


### Safety

During the pivotal OASIS trial ITT analysis, Revi demonstrated a favorable safety profile at 12 months.<sup>1</sup>




### Outcomes at 12 Months

in the Completers Analysis<sup>†</sup> (N=139) Include<sup>1</sup>:



### Outcomes at 12 Months

in the Intent-To-Treat (ITT) Analysis Include<sup>1</sup>:



#### Reference:

1. Heesakkers, J., et al. Pivotal Study of a Novel Wirelessly Powered, Patient Tailored Programmed, Tibial Neurostimulator for the Treatment of Patients with Overactive Bladder. Journal of Urology. 2023; 209(4S):e1178.

\*Due to battery depletion, implant migration, or lead fracture.

<sup>†</sup>Completers analysis represents the patients with available data at 12 months. Satisfaction results are from 130 patients.



## The Power of Proof<sup>1</sup>

In the OASIS Study Completers analysis<sup>†</sup> (N=139) at 12 months:



## How it Works

### Proven Process. Single Procedure.

Revi is a patient-centric therapy with a unique implant placement and system design. It provides an alternative treatment solution that may be better suited to some of your UUI patients.

### Convenient Therapy.

Revi is designed to be easy to use at home so it can fit into patients' lives and help them stay compliant with therapy.



**Minimally Invasive:** The battery-free miniature device is implanted in the ankle region near the posterior tibial nerve in a single predictable procedure under local anesthesia.



**Lightweight Wearable:** 7.8 oz wireless Wearable is worn around the ankle to activate the implant, only during therapy.



**Predictable Placement:** Subfascial placement allows for direct visualization of the posterior tibial nerve. Intraoperative testing and suture fixation allow for predictable control.



**Convenient Sessions:** At the patient's convenience for about 30-60 minutes daily, based on symptoms.



**Customized Therapy:** Optimized treatment based on patient response. Clinician Programmer, Hub, and myRevi App allow provider visibility into device usage and status to better support and optimize patient therapy.



**Customized Programming:** With preset treatment programs set for use at home.

#### Reference:

1. Heesakkers, J., et al. Pivotal Study of a Novel Wirelessly Powered, Patient Tailored Programmed, Tibial Neurostimulator for the Treatment of Patients with Overactive Bladder. *Journal of Urology*. 2023; 209(4S):e1178.

<sup>†</sup>Completers analysis represents the patients with available data at 12 months. Satisfaction results are from 130 patients.

# Designed With Both Patients and the Physicians Who Treat Them in Mind

Integrated components make the procedure, treatment, and ongoing management convenient.



## Revi Implant

- Advanced battery-free implant allows the device to be miniature in size without the need for replacement surgeries due to battery depletion, lead fracture, or lead migration.
- Only 3 cm in length and 3 mm in diameter, this tibial implant delivers reliable and long-lasting performance in a compact form.
- MR Conditional



## Revi Wearable

- Innovative lightweight external wearable technology powers the implant, enabling programming changes without regard for energy consumption.
- A mere 7.8 oz, it is worn around the ankle for ~30-60 minutes daily based on symptoms, delivering therapy at the patient's convenience.



## Revi Clinician Programmer

- Allows for customization of therapy based on a patient's response.
- Enables adjustment of frequency, pulse width, polarity, and amplitude.



## Hub

- Wirelessly retrieves therapy data from the Wearable.
- Provides the clinical care team visibility into device usage and status to better support and optimize patient therapy.



## myRevi App

- Online patient support center providing customized education, online chat functionality, and symptom tracking.
- Enables monitoring of symptoms to facilitate therapy optimization.

## Learn About Our Comprehensive Services For:



**Patient Support Program**



**Clinician Training**

### Indication

Revi is indicated for the treatment of patients with symptoms of urgency incontinence alone or in combination with urinary urgency.

Visit [www.bluewindmedical.com](http://www.bluewindmedical.com) for more details.

### Caution:

Federal Law (USA) restricts this device TO sale by or on the order of a physician. Refer to product instruction manual/package insert for instructions, warning, precautions, and contraindications. For more information, please contact BlueWind Medical at 888-715-2080, and/or consult BlueWind Medical's website at [www.bluewindmedical.com](http://www.bluewindmedical.com).

